



LABORATORY SERVICES

TEST CHANGES EFFECTIVE SEPTEMBER 13th, 2004

Please note the following test updates. Make the changes in your Interpath Service Manual if appropriate. Contact Client Services at 800-700-6891 for additional information.

<u>Test Number</u>	<u>Test Name</u>	<u>Change</u>
931042 (Also as part of Test number 91234 Von Willibrands Panel)	Factor VIII Activity	New reference range: 56-191%
91332 & 91327	GHB, Serum and Urine	New CPT code: 80100
91140	IGF BP3	New reporting units: ng/mL New reference range: Extensive age/sex ranges. Stability: Submit frozen, stable refrigerated 24-hours.

AGE	MALE	FEMALE
0-1 year	722-2204 ng/mL	588-2458 ng/mL
2-3 years	793-2222 ng/mL	935-2713 ng/mL
4-5 years	853-2541 ng/mL	1196-2468 ng/mL
6-7 years	1069-2234 ng/mL	1494-2498 ng/mL
8-9 years	1246-2900 ng/mL	1646-3102 ng/mL
10-11 years	1273-3777 ng/mL	1407-3579 ng/mL
12-13 years	1485-3668 ng/mL	1967-3501 ng/mL
14-15 years	1739-3417 ng/mL	1922-3421 ng/mL
16-17 years	1652-3133 ng/mL	1891-3202 ng/mL
18-19 years	1470-3316 ng/mL	1670-3246 ng/mL
20-24 years	1302-2974 ng/mL	1616-2996 ng/mL
25-29 years	1407-2896 ng/mL	1663-2929 ng/mL
30-34 years	1250-2903 ng/mL	1439-3325 ng/mL
35-39 years	1237-2604 ng/mL	1393-3042 ng/mL
40-44 years	1180-2780 ng/mL	1257-3007 ng/mL
45-49 years	1157-2850 ng/mL	1419-2477 ng/mL
50-54 years	1264-2525 ng/mL	1281-2798 ng/mL
55-59 years	1241-2730 ng/mL	1287-2957 ng/mL
60-64 years	1296-2385 ng/mL	1342-2565 ng/mL
65 years and older	1349-2814 ng/mL	1231-2637 ng/mL

91167	Mumps, IgG	New reference range: 0.89 IV or less = Negative 0.90-1.09 IV = Equivocal 1.10 IV or greater = Positive
-------	------------	---

Interpretation: IgG antibody to Mumps IgG detected, indicating a current or previous exposure/immunization. Positive IgG antibody levels in absence of current clinical symptoms may indicate immunity.

91187	Protein C, Functional	New reference range: 59-147%
91190	Protein S, Functional	New reference range: 66-143%

Continued/...

LABORATORY SERVICES Continued

TEST CHANGES EFFECTIVE SEPTEMBER 13th, 2004

Please note the following test updates. Make the changes in your Interpath Service Manual if appropriate. Contact Client Services at 800-700-6891 for additional information.

<u>Test Number</u>	<u>Test Name</u>	<u>Change</u>
93370 (Also as part of Test number 91118 Heavy Metals Panel, Blood)	Mercury, Blood	New reference range: 0-10 ug/L
91220	Tricyclic Screen	New methodology: EIA New interpretation: Results represents the total level of all detectable tricyclic antidepressants present. Tricyclics detectable with this assay include: amitriptyline, clomipramine, desipramine, doxepin, imipramine, nortriptyline, protriptyline and trimipramine. Levels of tricyclic antidepressants in the therapeutic range may not be detectable. False positives may occur with the following drugs: Seroquel, Triptal, Benadryl at toxic concentrations, Flexeril, Thioridazine, and Thorazine.